

<b>STN</b>	<b>Postup posudzovania expozície pracovníkov s aktívnymi implantovateľnými zdravotníckymi pomôckami elektromagnetickým poľom Časť 2-2: Špecifické hodnotenie pre pracovníkov s implantovateľnými kardioverter defibrilátormi (ICD)</b>	<b>STN EN 50527-2-2</b>  36 7938
------------	--	--

Procedure for the assessment of the exposure to electromagnetic fields of workers bearing active implantable medical devices - Part 2-2:  
Specific assessment for workers with cardioverter defibrillators (ICDs)

Táto norma obsahuje anglickú verziu európskej normy.  
This standard includes the English version of the European Standard.

Táto norma bola oznámená vo Vestníku ÚNMS SR č. 09/18

Obsahuje: EN 50527-2-2:2018

**127208**

EUROPEAN STANDARD

**EN 50527-2-2**

NORME EUROPÉENNE

EUROPÄISCHE NORM

May 2018

ICS 11.040.40; 17.240

English Version

Procedure for the assessment of the exposure to  
electromagnetic fields of workers bearing active implantable  
medical devices - Part 2-2: Specific assessment for workers with  
cardioverter defibrillators (ICDs)

Procédure pour l'évaluation de l'exposition des travailleurs  
porteurs de dispositifs médicaux implantables actifs aux  
champs électromagnétiques - Partie 2-2 : Evaluation  
spécifique aux travailleurs porteurs de défibrillateurs  
automatiques implantables

Verfahren zur Beurteilung der Exposition von  
Arbeitnehmern mit aktiven implantierbaren medizinischen  
Geräten gegenüber elektromagnetischen Feldern - Teil 2-2:  
Besondere Beurteilung für Arbeitnehmer mit Cardioverter-  
Defibrillatoren (ICDs)

This European Standard was approved by CENELEC on 2018-04-03. CENELEC members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration.

Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN-CENELEC Management Centre or to any CENELEC member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CENELEC member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CENELEC members are the national electrotechnical committees of Austria, Belgium, Bulgaria, Croatia, Cyprus, the Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, the Netherlands, Norway, Poland, Portugal, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom.



European Committee for Electrotechnical Standardization  
Comité Européen de Normalisation Electrotechnique  
Europäisches Komitee für Elektrotechnische Normung

**CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels**

**EN 50527-2-2:2018 (E)**

<b>Contents</b>	<b>Page</b>
European foreword .....	5
<b>1 Scope .....</b>	<b>6</b>
<b>2 Normative references .....</b>	<b>6</b>
<b>3 Terms and definitions .....</b>	<b>7</b>
<b>4 Specific assessment .....</b>	<b>8</b>
<b>4.1 Description of the assessment process .....</b>	<b>8</b>
4.1.1 General .....	8
4.1.2 Equipment consideration .....	11
4.1.3 Patient warning consideration .....	11
4.1.4 Cases for additional investigation .....	11
4.1.5 Choice of investigative method .....	14
4.2 Clinical investigation .....	15
4.3 Non-clinical investigation .....	16
4.3.1 General .....	16
4.3.2 Non-clinical investigation by <i>in vitro</i> testing .....	16
4.3.3 Non-clinical investigation by comparative study .....	19
<b>5 Documentation .....</b>	<b>20</b>
<b>Annex A (normative) Device specific replacement of EN 50527-1:2016, Table 1 .....</b>	<b>21</b>
<b>Annex B (informative) Clinical investigation methods .....</b>	<b>27</b>
B.1 Assessment of device compatibility using stored data and diagnostic features .....	27
B.2 Real time event monitoring by telemetry .....	27
<b>Annex C (informative) <i>in vitro</i> testing/measurements .....</b>	<b>28</b>
C.1 Introduction .....	28
C.2 EM phantom .....	28
C.2.1 General .....	28
C.2.2 EM phantom design .....	28
C.3 Basic procedure for device <i>in vitro</i> testing .....	29
<b>Annex D (informative) Modelling — Field modelling or calculations .....</b>	<b>31</b>
<b>Annex E (informative) Interference from Low-Frequency Magnetic and Electric Fields (1 Hz to 10 MHz) Including Application to 50 Hz Power-Frequency .....</b>	<b>32</b>
E.1 Introduction .....	32
E.2 Implanted devices and leads .....	32
E.3 Sensitivity of devices to interference .....	33
E.4 Immunity requirements .....	33
E.5 Voltage induced in the leads by magnetic fields .....	34
E.5.1 Induction in a loop .....	34
E.5.2 Inductive loop area for bipolar leads .....	34
E.5.3 Voltage induced in leads .....	36
E.6 Voltage induced in the leads by electric fields .....	36
E.7 Values of 50 Hz magnetic and electric field that can cause interference .....	38
E.8 Factors that affect the immunity from interference .....	38
E.8.1 Reasons for improved immunity .....	38
E.8.2 Adjustment for device sensitivity setting .....	39
E.8.3 Adjustment for lead tip to ring spacing .....	41
E.9 Application to Power Frequency Exposure Situations .....	42
E.9.1 Public exposures .....	42
E.9.2 Beneath high voltage power lines .....	42
E.9.3 Occupational settings .....	43
E.9.4 Temporary exposure above the interference levels .....	43
E.9.5 Induced voltages at occupational exposure limit values .....	44

<b>E.10</b>	<b>Conversion based on known compliance with basic restrictions .....</b>	<b>44</b>
E.10.1	General.....	44
E.10.2	Relationship between magnetic fields and induced current density.....	45
E.10.3	Relationship between magnetic fields and induced voltages on an implanted lead .....	46
E.10.4	A simple model to analyse the possible voltages at device terminations generated from induced current density equivalent the basic restrictions of Council Recommendation 1999/519/EC .....	46
<b>Annex F</b> (informative)	<b>Determination of minimum immunity for radio-frequency fields .....</b>	<b>49</b>
<b>F.1</b>	<b>Determination of immunity from the field .....</b>	49
F.1.1	General.....	49
F.1.2	Intermediate frequencies (5 MHz to 30 MHz) .....	49
F.1.3	High frequency (above 30 MHz) .....	50
<b>F.2</b>	<b>References .....</b>	51
<b>Annex G</b> (informative)	<b>Obtaining the device immunity and guidelines provided by device manufacturers – .....</b>	<b>52</b>
<b>G.1</b>	<b>Introduction.....</b>	52
<b>G.2</b>	<b>EMC and devices – General guidelines, hazards and harms.....</b>	52
<b>G.3</b>	<b>Induced voltages, fields and zones .....</b>	55
G.3.1	Induced voltage test levels .....	55
G.3.2	Magnetic field amplitudes producing test limits .....	55
G.3.3	Induced voltage zones .....	57
G.3.4	Magnetic field zones.....	58
Bibliography .....		60

## Figures

<b>Figure 1</b>	<b>— Overview of the assessment process.....</b>	<b>9</b>
<b>Figure 2</b>	<b>— Specific assessment process.....</b>	<b>10</b>
<b>Figure 3</b>	<b>— Additional investigation process .....</b>	<b>13</b>
<b>Figure 4</b>	<b>— Comparison process .....</b>	<b>18</b>
<b>Figure C.1</b>	<b>— Example of <i>in vitro</i> procedure for EM interference at low frequency using planar electrodes, bipolar lead and ECG and data recording.....</b>	<b>30</b>
<b>Figure E.1</b>	<b>— Typical implantations of cardiac devices .....</b>	<b>32</b>
<b>Figure E.2</b>	<b>— Effective induction area of an open wire loop inside a conductive medium .....</b>	<b>35</b>
<b>Figure E.3</b>	<b>— Schematic representation of bipolar pickup of interference in an infinitely homogeneous conducting medium.....</b>	<b>35</b>
<b>Figure E.4</b>	<b>— Showing how the immunity ratio affects magnetic field that can result in interference .....</b>	<b>40</b>
<b>Figure E.5</b>	<b>— Showing how the immunity ratio affects electric field that can result in interference ...</b>	<b>40</b>
<b>Figure E.6</b>	<b>— Showing how the tip to ring spacing affects the magnetic field that can result in interference .....</b>	<b>41</b>
<b>Figure E.7</b>	<b>— Showing how the tip to ring spacing affects the electric field that can result in interference .....</b>	<b>42</b>
<b>Figure E.8</b>	<b>— Eddy-current inside a conductive medium induced by varying magnetic flux .....</b>	<b>45</b>
<b>Figure E.9</b>	<b>— Voltage induced on a lead inside conductive body tissue .....</b>	<b>46</b>
<b>Figure E.10</b>	<b>— Voltages on an implanted lead.....</b>	<b>47</b>
<b>Figure G.1</b>	<b>— Induced voltage test levels .....</b>	<b>55</b>
<b>Figure G.2</b>	<b>— Magnetic field amplitudes, for frequencies below 10 000 kHz, producing test limits in bipolar configurations.....</b>	<b>56</b>
<b>Figure G.3</b>	<b>— Induced voltage zones for bipolar configurations.....</b>	<b>57</b>
<b>Figure G.4</b>	<b>— Induced voltage zones for bipolar configurations.....</b>	<b>58</b>

**EN 50527-2-2:2018 (E)**

**Figure G.5 — Magnetic field zones, for frequencies below 10 000 kHz for bipolar configurations ..... 59**

**Tables**

<b>Table A.1 — Compliant workplaces and equipment with exceptions .....</b>	<b>21</b>
<b>Table E.1 — Values of 50 Hz electric and magnetic field (r.m.s.) that might, under unfavourable circumstances, cause interference with devices .....</b>	<b>38</b>
<b>Table E.2 — Summary of typical maximum field values beneath high-voltage overhead lines and corresponding voltage induced in the device .....</b>	<b>43</b>
<b>Table E.3 a) — Voltage induced in bipolar leads by the health effects ELV where f is in Hz .....</b>	<b>44</b>
<b>Table E.3 b) — Voltage induced in bipolar leads by the sensory effects ELV, where f is in Hz. It is defined only up to 400 Hz.....</b>	<b>44</b>

## European foreword

This document (EN 50527-2-2:2018) has been prepared by CLC/TC 106X "Electromagnetic fields in the human environment".

The following dates are fixed:

- latest date by which this document has to be implemented at national level by publication of an identical national standard or by endorsement (dop) 2019-04-03
- latest date by which the national standards conflicting with this document have to be withdrawn (dow) 2021-04-03

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CENELEC shall not be held responsible for identifying any or all such patent rights.

This document has been prepared under a mandate given to CENELEC by the European Commission and the European Free Trade Association.

**EN 50527-2-2:2018 (E)****1 Scope**

This European Standard provides the procedure for the specific assessment required in EN 50527-1:2016, Annex A, for workers with implanted cardioverter defibrillators (ICDs) and Cardiac Resynchronization Therapy devices with associated defibrillation functions (CRT-D). Only devices of this type equipped with leads implanted transvenously are considered. It offers different approaches for doing the risk assessment.

NOTE 1 If the worker has other Active Implantable Medical Devices (AIMDs) implanted additionally, they are assessed separately according to EN 50527-1 or other particular standards within the EN 50527 series.

NOTE 2 The risks to patients due to interference with pacing functions associated with CRT-D devices are assessed using EN 50527-2-1.

The purpose of the specific assessment is to determine the risk for workers with implanted ICDs and CRT-Ds arising from exposure to electromagnetic fields (EMF) at the workplace. The assessment includes the likelihood of clinically significant effects and takes account of both transient and long-term exposure within specific areas of the workplace.

NOTE 3 This standard does not address risks from contact currents.

The techniques described in the different approaches may also be used for the assessment of publicly accessible areas.

The frequency range to be observed is from 0 Hz to 3 GHz. Above 3 GHz no interference with the devices within the scope of this Particular Standard is expected to occur when the exposure limits are not exceeded.

NOTE 4 The rationale for limiting the observation range to 3 GHz can be found in ISO 14117:2012, Clause 5.

NOTE 5 Further information concerning the functions of Pacemakers, CRT-D, and ICD devices can be found in Ellenbogen and Kaszala, 2014.

**2 Normative references**

The following documents, in whole or in part, are normatively referenced in this document and are indispensable for its application. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN 45502-2-2:2008, *Active implantable medical devices — Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)*

EN 50527-1:2016, *Procedure for the assessment of the exposure to electromagnetic fields of workers bearing active implantable medical devices — Part 1: General*

EN ISO 14155, *Clinical investigation of medical devices for human subjects — Good clinical practice (ISO 14155)*

ISO 14117:2012, *Active implantable medical devices — Electromagnetic compatibility — EMC test protocols for implantable cardiac pacemakers, implantable cardioverter defibrillators and cardiac resynchronization devices*

**koniec náhľadu – text ďalej pokračuje v platenej verzii STN**